

coldes

Textilien GmbH


Sterile Surgical Gown

FABRIC
SMS 40 GSM NONWOVEN
Disposable

Color : Blue - White
Standart : CE 2163 / EN 13795 / EN 13485
Category : Level 2

USAGE

Used in **LOW** risk situations



PENETRATION

Provides a barrier to larger amounts of fluid penetration through splatter and some fluid exposure through soaking two tests are conducted to test barrier protection performance:

- 1 - Water impacting the surface of the gown material
- 2- Pressurizing the material

APPLICABILITY

Basic care, Standart hospital medical unit



Weight
40 gr/m²



CE
2163

EN 13795

ISO
13485

INTPE02-S